

# Liposuction Cannula Catalogue

Manufacturers & Exporters



**TETRA INDUSTRIES**

SURGICAL & MEDICAL INSTRUMENTS

## Liposuction Cannula



# Liposuction Cannula Catalogue

## Liposuction Cannula Handles



### BLACK HANDLE

Product Code:AK-LP-10 Product Description:

Black Handle

Handle made of aluminum with type 2 black anodized, 6.65" length 1" diameter



### BLUE HANDLE

AK-LP-15

Product Description:

BLUE Handle

Handle made of aluminum with type 2 blue anodized, 6.65" length 1" diameter



### GOLD HANDLE

Product Code:AK-LP-20

Product Description:

Gold Handle

Handle made of aluminum with type 2 gold anodized, 6.65" length 1" diameter



### HANDLE w/SWIVEL BARB

Product CodeA:K-LP-25

Product Code:

Product Description:

Handle made of aluminum with type 2 blue anodized, 6.65" length, 1" diameter for Large Handle, 3/4" for Small Handle)



### HANDLE

Product Code:AK-LP-30 Product Code:AK-LP-35

Product Description:

Handle made of aluminum with type 2 blue anodized, 6.65" length 1" diameter for Large Handle, 3/4" for Small Handle



TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Liposuction Cannula Handles Cont.



### STAINLESS STEEL TOOMEY 60cc HUB

Product Code:AK-LP-50

**Product Description:**

Stainless Steel and Stainless Steel hub



### ALUMINUM TOOMEY 60cc HUB

Product CodeA:K-LP-55

**Product Description:**

Hub made of Aluminum with type 2 blue anodized, 1.55" length 1.31" diameter



### QUICK SCREW HUB

Product CodeA:K-LP-60

**Product Description:**

Hub made of Stainless Steel, 2.1" length 7/8" diameter



### SMALL QUICK DISCONNECT HUB

Product CodeA:K-LP-65

**Product Description:**

Hub made of Stainless Steel, 1.93" length 5/8" diameter



### QUICK DISCONNECT HUB

Product CodeA:K-LP-70

**Product Description:**

Hub made of Stainless Steel, 1.6" length 7/8" diameter



### QUICK LOCK HUB

Product CodeA:K-LP-75

Product Description: Hub made of Stainless Steel, 1.65" length 1/2" diameter



### LUER LOCK HUB

Product CodeA:K-LP-80

**Product Description:**

Luer Lock Hub

Hub made of Stainless Steel, 1.16" length 3/8" diameter



TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Cannula Hole Pattern & Cont.



### MERCEDES

Product Code: AK-LP-80  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### TRI PORT

Product Code: AK-LP-85 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### LAS VEGAS

Product Code: AK-LP-90 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### STANDARD 3 PORT

Product Code: AK-LP-95 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### STANDARD 2 PORT

Product Code: AK-LP-100 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm. Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### STANDARD SINGLE PORT

Product Code: AK-LP-105 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### 3 PORT 60 DEGREE INDEX

Product Code: AK-LP-110 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### 3 PORT 60 DEGREE INDEX

Product Code: AK-LP-115 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### BASKET

Product Code: AK-LP-120 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### SPATULA

Product Code: AK-LP-125 Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.

TETRA INDUSTRIES

# Liposuction Cannula Catalogue

## Cannula Hole Pattern & Cont.



### SPATULA 2 PORT

Product Code:AK-LP-130  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### SPATULA 3 PORT

Product Code:AK-LP-135  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### SPATULA 3 PORT

Product Code:AK-LP-140  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### COBRA

Product Code:AK-LP-145  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### FUENTES

Product Code:AK-LP-150  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### SAYLAN

Product Code:AK-LP-155  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### CANDY CANE

Product Code:AK-LP-160  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### MODERN SPIRAL

Product Code:AK-LP-165  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### TOLEDO SINGLE PORT

Product Code:AK-LP-170  
Available Diameter & Length  
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Infiltration & Reinjection Micro Cannula



### GARDEN SPRAY ANESTHESIA NEEDLE

Product Code:AK-LP-175  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga,  
3mm, 4mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### MULTI HOLE SPRAY NEEDLE

Product Code:AK-LP-180  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga,  
3mm, 4mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### BODY ANESTHESIA NEEDLE

Product Code:AK-LP-185  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga,  
3mm, 4mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### Infusion 6 hole Needle

Product Code:AK-LP-190  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga,  
3mm, 4mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### Infusion 8 hole Needle

Product Code:AK-LP-195  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga,  
3mm, 4mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



### MICRO INJECTION NEEDLE, STYLE 3

Product Code: AK-LP-200  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm,  
4mm.  
Length: 15cm, 23cm, 26cm.



### MICRO INJECTION NEEDLE, STYLE 2

Product Code: AK-LP-205  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm,  
4mm.  
Length: 15cm, 23cm, 26cm.



### MICRO INJECTION NEEDLE, STYLE 1

Product Code: AK-LP-210  
Available Diameter & Length  
Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm,  
4mm.  
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



# Liposuction Cannula Catalogue

## Reinjection Micro Cannula & Cont.



### MICRO INJECTION NEEDLE "V" DISSECTOR

Product Code:AK-LP-215

Available Diameter & Length

Diameter:18ga,17ga,16 ga.  
Length: 5cm, 7cm, 8cm, 9cm, 12cm.



### BEVEL 45

Product Code:AK-LP-220

Available Diameter

Length Diameter: 16ga,  
15ga, 14ga.



### BEVEL 30

Product Code:AK-LP-225

Available Diameter

Length Diameter: 16ga,  
15ga, 14ga.



**TETRA INDUSTRIES**  
SURGICAL & MEDICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Infusion Handles & (Cont.)



### INFUSION HANDLE WITH ON/OFF CONTROL

Model Number: AK-LP-500

Product Description:

A flow on/off control infusion handle for use with roller pump devices.

Handle diameter: 15/16" (Stainless Steel control)

Valve: Stainless Steel

Female Luer Lock: Stainless Steel

Male Luer Lock: Stainless Steel



### INFUSION HANDLE WITH ON/OFF CONTROL

Model Number: AK-LP-550

Product Description:

A flow on/off control infusion handle for use with roller pump devices.

Handle diameter: 15/16" (control) 3/4" (handle)

Valve: Stainless Steel

Female Luer Lock: Stainless Steel

Male Luer Lock: Stainless Steel



### INFUSION HANDLE WITH ON/OFF STOPCOCK CONTROL

Model Number: AK-LP-600

Product Description:

Infusion Handle with Stopcock Control

Handle diameter: 3/4"

Valve: Brass with Electroless Nickel Plating

Female Luer Lock: Stainless Steel

Male Luer Lock: Stainless Steel



### INFUSION HANDLE WITH ON/OFF STOPCOCK CONTROL

Model Number: AK-LP-650

Product Description:

Infusion Handle with Stopcock Control

Handle diameter: 3/4"

Valve: Brass with Electroless Nickel Plating

Female Luer Lock: Stainless Steel

Male Luer Lock: Stainless Steel



### INFUSION HANDLE WITH SWIVEL LUER NUT

Model Number:

AK-LP-700

AK-LP-750 (for 1" Dia.)

Product Description:

Standard Infusion Handle, No Control with Swivel Nut

Handle diameter: 3/4"

Female Luer Lock: Stainless Steel

Male Swivel Luer Lock: Stainless Steel



### LUER LOCK ADJUSTABLE HANDLE

Model Number:

AK-LP-800 (for small)

AK-LP-850 (for large)

Product Description:

Luer Lock Adjustable Handle. Orientation between Luer Lock and Handle.

Handle diameter: 3/4" (for small) 1.0" (for large)

Female Luer Lock: Stainless Steel

Male Luer Lock: Stainless Steel



# Liposuction Cannula Catalogue

## Infusion Handles & (Cont.)



### INFUSION HANDLE WITH TUBE CONNECTOR

Model Number:  
AK-LP-900  
AK-LP-950 (for 1" Dia.)

Product Description:  
Infusion handle with tube connector for use with reusable infusion tubing.  
Handle diameter: 3/4"  
Female Luer Lock: Stainless Steel  
Barb: Stainless Steel



### LUER LOCK ADJUSTABLE ASPIRATION HANDLE

Model Number: AK-LP-1000

Product Description:  
Luer Lock Adjustable Handle ideal for rotating the Luer Lock flat  
Handle diameter: 3/4"  
Male Luer Lock: Stainless Steel



### ASPIRATION HANDLE

Model Number:  
AK-LP-1050 (for large)  
AK-LP-1100 (for small)

Product Description:  
Standard Aspiration Handle  
Handle diameter: 3/4" (for small) 1.0" (for large)  
Male Luer Lock: Stainless Steel



### INFUSION HANDLE WITH SWIVEL LUER LOCK

Model Number:  
AK-LP-1150  
AK-LP-1200L (for 1" Dia.)

Product Description:  
Standard Infusion Handle, No Control with Swivel Nut  
Handle diameter: 3/4"  
Female Luer Lock: Stainless Steel  
Male Swivel Luer Lock: Stainless Steel

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Two Piece Cannuals & (Cont.)



### QUICK DISCONNECT CANNULA SET

Model Number:AK-LP-1250

Product Description:

A quick Disconnect Handle is a quick and easy interchangeable cannula hub.

Overall length: 6.65"

Body: All Stainless Steel



### QUICK LOCK CANNULA SET

Model Number:AK-LP-1300

Product Description:

A quick Lock Handle is a quick and easy interchangeable Cannulas hub.

Standard vented holes and swivel barb

Overall length: 6.65"

Body: Aluminum with type 2 blue anodized

Lock : Stainless Steel

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS



### QUICK LOCK CANNULA SET

Model Number:AK-LP-1350

Product Description:

A quick Lock Handle is a quick and easy interchangeable Cannulas hub.

Standard vented holes and swivel barb

Overall length: 6.65"

Body: Aluminum with type2 blue anodized

Lock : Stainless Steel



### TOOMEY ASPIRATION CANNULA SET

Model Number:AK-LP-1400

Product Description:

Toomey Aspiration Handle

Overall length: 7.1"

Body: Delrin (Autoclavable)

# Liposuction Cannula Catalogue

## Two Piece Cannulas & (Cont.)

### LUER ADJUSTABLE

Part Number: AHKA-LNP-D14L5E0  
LOCK



Product Description:  
Luer Lock Adjustable Handle  
ideal for rotating the Luer Lock  
flat to the thumb rest  
Handle diameter: 3/4"  
Overall length: 7.2"  
Body: Aluminum with type 2  
blue anodized  
Male Luer Lock: Stainless Steel

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS



### ASPIRATION HANDLE

Part Number:  
AK-LP-1500 (for large)  
AK-LP-1550 (for small)

Product Description:  
Standard Aspiration Handle  
Handle diameter: 3/4" (for small)  
1.0" (for large)  
Overall length: 6.85"  
Body: Aluminum with type 2 blue  
anodized  
Male Luer Lock: Stainless Steel

# Liposuction Cannula Catalogue

## Injection Gun Set



### INJECTION GUN SET 1ML, 3ML, 5ML, 10ML BD

Model Number:AK-LP-2000

#### Product Description:

Standard method of micro tissue injections. All stainless steel

For 1ml, 3ml, 5ml, 10ml BD syringes  
.15 ml per trigger on 1ml BD syringe  
.32 ml per trigger on 3ml BD syringe  
.62 ml per trigger on 5ml BD syringe  
.88 ml per trigger on 10ml BD syringe



### INJECTION GUN SET 3CC, 5CC, 10CC BD

Model Number:AK-LP-3000 Product Description:

Standard method of micro tissue injection. All stainless steel

For 3cc, 5cc, 10cc BD syringes

.32 ml per trigger on 3ml BD syringe  
.62 ml per trigger on 5ml BD syringe  
.88 ml per trigger on 10ml BD syringe

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS



### INJECTION GUN SET 3CC

Model Number:AK-LP-4000

#### Product Description:

Standard method of micro tissue injection. All stainless steel

For 3cc BD syringes

.32 ml per trigger on 3ml BD syringe



### INJECTION GUN SET 5CC

Model Number:AK-LP-5000

#### Product Description:

Standard method of micro tissue injection. All stainless steel

For 5cc BD syringes

.62 ml per trigger on 5ml BD syringe

# Liposuction Cannula Catalogue

## Injection Gun Set



### INJECTION GUN SET 10CC

Model Number:AK-LP-6000

Product Description:

Standard method of micro tissue injection. All stainless steel  
For 10cc BD syringes  
.88 ml per trigger on 10ml BD syringe



### INJECTION GUN SET 20CC

Model Number: AK-LP-7000

Product Description:

Standard method of micro tissue injection. All stainless steel  
For 20cc BD syringes  
1.5 ml per trigger on 20ml BD syringe



### INJECTION GUN SET 30CC

Model Number:AK-LP-8000

Product Description:

Standard method of micro tissue injections. All stainless steel  
For 30cc BD syringes  
1.7 ml per trigger on 30ml BD syringe



### INJECTION GUN SET 60CC

Model Number:AK-LP-9000

Product Description:

Standard method of micro tissue injections. All stainless steel  
For 60cc Monoject syringes  
1.9 ml per trigger on 60ml Monoject syringe

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Injection Gun Set



### Additional Syringe Accessories

Part Number:	Product Description:
AK-LP-	Plunger Adaptor Removal
1500 AK-	Plunger adaptor for 3 ml
LP-1550	Plunger adaptor for 5 ml
AK-LP-	Plunger adaptor for 10 ml
1560 AK-	Plunger adaptor for 20 ml
LP-1570	Plunger adaptor for 30 ml
AK-LP-	Plunger adaptor for 35 ml
1580 AK-	Plunger adaptor for 60 ml

LP-1590

AK-LP-

1600 AK-

LP-1650

Product Description:

This additional accessories are ideal for a surgeon who wants to have two or more syringes ready for tissue re injection procedure.

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Toomey & Accessories



### TOOMEY TO LUER LOCK SYRINGE

Model Number:AK-LP-2002

Product Description:  
Toomey 60cc Monoject to Luer Lock Syringe  
Stainless Steel and  
Stainless Steel hub



### TOOMEY TO LUER LOCK SYRINGE

Model Number:AK-LP-2003

Product Description:  
Toomey 60cc Monoject to Luer Lock Syringe  
Female Luer Lock made of Stainless Steel and  
Aluminum hub



### TOOMEY TO LUER LOCK NEEDLE

Model Number:AK-LP-2004

Product Description:  
Toomey 60cc Monoject to Luer Lock  
Needle Stainless Steel  
hub

TETRA INDUSTRIES  
SURGICAL INSTRUMENTS

# Liposuction Cannula Catalogue

## Toomey & Accessories



### TOOMEY TO LUER LOCK NEEDLE

Model Number:AK-LP-1450

Product Description:  
Toomey 60cc Monoject to Luer Lock Needle  
Male Luer Lock made of Stainless Steel and  
Aluminum hub



### SYRINGE BARB ADAPTOR

Part Number:  
AK-LP-1460 SBA60SS  
AK-LP-1465 SBA35SS  
AK-LP-1462 SBA10SS

Product Description:  
Syringe Barb Adapter  
Made of Stainless Steel



### LUER LOCK TO LUER LOCK

Model Number:AK-LP-1468

Product Description:  
Luer Lock syringe to Luer Lock syringe soft Tissues  
Transfer system  
Made of Stainless Steel

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS



### Sterilization Tray



#### STERILIZATION TRAY

Model No AK-LP-3600 7trays are sterilizable by all standard methods of

sterilization. Trays are designed to encase instruments during processing. Safe for all standard sterilization methods. Autoclave/dry heat up to 360°F, ETO, and Cold Solutions.

Flexibility of layout. Fingertip mats.

Unique, "grid" system in all bases, trays, and lids makes it easy to install mats or inserts to protect your instruments.

4.0x7.5x1.5 Micro Tray Base Lid & Mat

TETRA INDUSTRIES  
SURGICAL & MEDICAL INSTRUMENTS



**TETRA INDUSTRIES**  
SURGICAL & MEDICAL INSTRUMENTS

**PHONE: +92-3297598399**

**Email:**

**tetraindustries01@gmail.com**